

Adipokinetic Hormone Peptide

Catalogue No.:BTA11171

Adipokinetic Hormone is a peptide.

Target:	Adipokinetic Hormone
Conjugation:	Unconjugated
Form:	Lyophilized
Purity:	> 95% (HPLC)
Storage:	Store dry at -20 °C. Avoid repeated freeze/thaw cycles.
Molecular Weight:	921 Da
Molecular Formula:	$C_{44}H_{60}N_{10}O_{12}$
Sequence:	Glp-Leu-Thr-Phe-Thr-Ser-Ser-Trp-Gly-NH ₂
Note:	This product is for research use only. Not for human consumption, cosmetic, therapeutic or diagnostic use.