

## Alsevalimab

Catalogue No.:BTA14520

This product is for research use only (RUO) by laboratories. It is not suitable for use in humans or as a medicine. Alsevalimab is a Human Monoclonal antibody targeting V-set domain-containing T-cell activation inhibitor 1, expressed in Mammalian cells.

<b>Target:</b>	V-set domain-containing T-cell activation inhibitor 1
<b>Clonality:</b>	Monoclonal
<b>Reactivity:</b>	Human
<b>Host:</b>	CHO cells
<b>Isotype:</b>	IgG <sub>1</sub> Kappa
<b>Form:</b>	Liquid
<b>Purity:</b>	> 95%
<b>Purification:</b>	Protein A or Protein G affinity purified.
<b>Storage:</b>	Store at 4°C for up to 2 weeks. Aliquot and store at -20°C or -80°C for up to 12 months.
<b>UniProt Primary AC:</b>	Q7Z7D3 ( <a href="#">UniProt</a> , <a href="#">Expasy</a> )
<b>Buffer:</b>	PBS, pH 7.4
<b>Concentration:</b>	1 mg/ml
<b>Note:</b>	This product is for research use only. Not for human consumption, cosmetic, therapeutic or diagnostic use.