

## BPP 9a Peptide

Catalogue No.:BTA11172

BPP 9a is a peptide.

<b>Target:</b>	BPP 9a
<b>Conjugation:</b>	Unconjugated
<b>Form:</b>	Lyophilized
<b>Purity:</b>	> 95% (HPLC)
<b>Storage:</b>	Store dry at -20 °C. Avoid repeated freeze/thaw cycles.
<b>Molecular Weight:</b>	1101.3 Da
<b>Molecular Formula:</b>	$C_{53}H_{76}N_{14}O_{12}$
<b>Sequence:</b>	Pyr-Trp-Pro-Arg-Pro-Gln-Ile-Pro-Pro
<b>Note:</b>	This product is for research use only. Not for human consumption, cosmetic, therapeutic or diagnostic use.