

Brucella Antibody (Anti-Brucella) Rapid Test Kit

Catalogue No.:BTA11112

Brucella Antibody Rapid Test Kit is a qualitative colloidal gold-based lateral flow immunochromatographic assay for the detection of Brucella Antibody.

Target: Brucella Antibody (Anti-Brucella)

Reactivity: Cow, Goat, Sheep, Pig

Tested Applications:LF

Storage: Store at 2-30 °C.

Validity: The validity for this kit is 24 months.

Detection Method: Immunochromatographic

Assay Type: Indirect

Assay Data: Qualitative

Sample Type: Serum and milk.

Note: This product is for research use only.

Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.