

CD3 (FITC) / CD4 (PE) Combination Antibody

Catalogue No.:BTA11201

CD3 (FITC) / CD4 (PE) Combination Antibody is an antibody cocktail consisting of two monoclonal antibodies:

- CD3 Mouse Monoclonal Antibody conjugated to FITC
- CD4 Rat Monoclonal Antibody conjugated to PE

Target:	CD3/CD4
Clonality:	Monoclonal, Monoclonal
Reactivity:	Dog
Tested Applications:	FCM
Host:	Mouse, Rat
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Conjugation:	FITC, PE
Isotype:	IgG ₁ , IgG _{2a}
Form:	Lyophilized
Storage:	Store at -20 °C or -80 °C. Avoid repeated freeze/thaw cycles.
UniProt Primary AC:	P27597 (UniProt , ExPAS) P33705 _____, ExPAS) y (UniProt) y
Gene Symbol:	CD3E, CD4
Buffer:	PBS, 1% BSA. Contains sodium azide.
Note:	This product is for research use only.