

CD4 / CD8 / CD3 Antibody (FITC, PE, PE / Cyanine 5)

Catalogue No.:BTA11204

CD4/CD8/CD3 Antibody (FITC/PE/PE-Cyanine 5) is a Mouse Monoclonal combination antibody against CD4/CD8/CD3.

Target:	CD4/CD8/CD3
Clonality:	Monoclonal
Reactivity:	Human
Tested Applications:	FCM
Host:	Mouse
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Conjugation:	FITC, PE, PE / Cyanine 5
Isotype:	IgG _{2b} , IgG ₁ , IgG _{2a}
Form:	Liquid
Storage:	Aliquot and store at -20°C. Avoid repeated freeze/thaw cycles.
UniProt Primary AC:	P01730 (UniProt , Expasy)
KEGG:	hsa:920
String:	9606.ENSP00000011653
Buffer:	PBS, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA.
Note:	This product is for research use only.