

CD4 (FITC) / CD8 (PE) Combination Antibody

Catalogue No.:BTA11203

CD4 (FITC) / CD8 (PE) Combination Antibody is an antibody cocktail consisting of two monoclonal antibodies:

- CD4 Antibody conjugated to FITC
- CD8 Antibody conjugated to PE

Target:	CD4/CD8
Clonality:	Monoclonal, Monoclonal
Reactivity:	Dog
Tested Applications:	FCM
Host:	Rat
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Conjugation:	FITC, PE
Isotype:	IgG _{2a} , IgG ₁
Form:	Lyophilized
Storage:	Store at -20 °C or -80 °C. Avoid repeated freeze/thaw cycles.
UniProt Primary AC:	P33705 (UniProt , ExPAS) P33706 _____, ExPAS) y (UniProt y)
Gene Symbol:	CD4, CD8A
Buffer:	PBS, 1% BSA. Contains sodium azide.
Note:	This product is for research use only.