

CD4 (FITC) / CD8 (PE) Combination Antibody

Catalogue No.:BTA11200

CD4 (FITC) / CD8 (PE) Combination Antibody is an antibody cocktail consisting of two monoclonal antibodies:

- CD4 Antibody conjugated to FITC
- CD8 Antibody conjugated to PE

Target: CD4/CD8

Clonality: Monoclonal, Monoclonal

Reactivity: Mouse

Tested Applications: FCM

Host: Rat

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Conjugation: FITC, PE

Isotype: IgG_{2b}, IgG_{2a}

Form: Lyophilized

Storage: Store at -20 °C or -80 °C. Avoid repeated freeze/thaw cycles.

UniProt Primary AC: P01731 ([UniProt](#), [ExPAS](#)) P06332 _____, [ExPAS](#)) P10300 _____, [ExPAS](#))
y (UniProt y (UniProt y

Gene Symbol: CD4, CD8A, CD8B1

Buffer: PBS, 1% BSA. Contains sodium azide.

Note: This product is for research use only.