

## Ciprofloxacin (CPFX) Rapid Test Kit

Catalogue No.:BTA11103

Ciprofloxacin (CPFX) Rapid Test Kit is a qualitative colloidal gold-based lateral flow immunochromatographic assay for the detection of Ciprofloxacin (CPFX).

**Target:** Ciprofloxacin (CPFX)

**Reactivity:** General

**Tested Applications:**LF

**Storage:** Store at 2°C to 30°C upon receipt.

**Validity:** The validity for this kit is 12 months.

**Sensitivity:** 1 ppb

**Detection Method:** Immunochromatographic

**Assay Type:** Competitive

**Assay Data:** Qualitative

**Sample Type:** Tissue, honey, egg and milk.

**Note:** This product is for research use only.

Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.