

Datasheet

Revision date: 08 Oct
2024

Crefmirlimab

Catalogue No.:BTA14505

This product is for research use only (RUO) by laboratories. It is not suitable for use in humans or as a medicine. Crefmirlimab is a Chimeric Humanized Monoclonal antibody, expressed in Mammalian cells.

Target: Crefmirlimab

Clonality: Monoclonal

Reactivity: Human

Host: CHO cells

Form: Liquid

Purity: > 95%

Purification: Protein A or Protein G affinity purified.

Storage: Store at 4°C for up to 2 weeks. Aliquot and store at -20°C or -80 °C for up to 12 months.

Buffer: PBS, pH 7.4

Concentration: 1 mg/ml

Note: This product is for research use only.
Not for human consumption, cosmetic, therapeutic or diagnostic use.