

Fibroblast Growth Factor 1, Acidic / AFGF (FGF1) Antibody (Biotin)

Catalogue No.:BTA11250

FGF1 Antibody (Biotin) is a Rabbit Polyclonal antibody against FGF1 conjugated to Biotin.

Target:	Fibroblast Growth Factor 1, Acidic / AFGF (FGF1)
Clonality:	Polyclonal
Reactivity:	Human
Tested Applications:	ELISA
Host:	Rabbit
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Conjugation:	Biotin
Immunogen:	Recombinant human Fibroblast growth factor 1 protein (16-155AA).
Isotype:	IgG
Form:	Liquid
Purity:	> 95%
Purification:	Purified by Protein G.
Storage:	Aliquot and store at -20°C. Avoid exposure to light. Avoid repeated freeze/thaw cycles.
UniProt Primary AC:	P05230 (UniProt , Expasy)
Gene Symbol:	FGF1
GeneID:	2246
OMIM:	131220
NCBI Accession:	NP_000791.1, NM_000800.4
KEGG:	hsa:2246

Datasheet

Revision date: 08 Oct 2024



BIOTECH ASSAY Inc.

String: [9606.ENSF00000480791](#)

Buffer: 0.01 M PBS, pH 7.4, 0.03% Proclin-300 and 50% Glycerol.

Note: This product is for research use only.

BIOTECH ASSAY