

## Fibroblast Growth Factor 7 (FGF7) Antibody (Biotin)

Catalogue No.:BTA11252

Fibroblast growth factor 7 Antibody (Biotin) is a Rabbit Polyclonal antibody against Fibroblast growth factor 7 conjugated to Biotin.

<b>Target:</b>	Fibroblast Growth Factor 7 (FGF7)
<b>Clonality:</b>	Polyclonal
<b>Reactivity:</b>	Human
<b>Tested Applications:</b>	ELISA
<b>Host:</b>	Rabbit
<b>Recommended dilutions:</b>	Optimal dilutions/concentrations should be determined by the end user.
<b>Conjugation:</b>	Biotin
<b>Immunogen:</b>	Recombinant human Fibroblast growth factor 7 protein (34-189AA).
<b>Isotype:</b>	IgG
<b>Form:</b>	Liquid
<b>Purity:</b>	> 95%
<b>Purification:</b>	Purified by Protein G.
<b>Storage:</b>	Aliquot and store at -20°C. Avoid exposure to light. Avoid repeated freeze/thaw cycles.
<b>UniProt Primary AC:</b>	P21781 ( <a href="#">UniProt</a> , <a href="#">Expasy</a> )
<b>Gene Symbol:</b>	FGF7
<b>GeneID:</b>	<a href="#">2252</a>
<b>OMIM:</b>	<a href="#">148180</a>
<b>NCBI Accession:</b>	NP_002000.1, NM_002009.3
<b>KEGG:</b>	hsa:2252

# Datasheet

Revision date: 08 Oct 2024



**BIOTECH ASSAY Inc.**

**String:** [9606.ENSP00000267843](#)

**Buffer:** 0.01 M PBS, pH 7.4, 0.03% Proclin-300 and 50% Glycerol.

**Note:** This product is for research use only.

BIOTECH ASSAY