

## Fin bud initiation factor homolog Antibody (Biotin)

Catalogue No.:BTA11257

Fin bud initiation factor homolog Antibody (Biotin) is a Rabbit Polyclonal antibody against Fin bud initiation factor homolog conjugated to Biotin.

**Target:** Fin bud initiation factor homolog

**Clonality:** Polyclonal

**Reactivity:** Human

**Tested Applications:** ELISA

**Host:** Rabbit

**Recommended dilutions:** Optimal dilutions/concentrations should be determined by the end user.

**Conjugation:** Biotin

**Immunogen:** Recombinant human Fin bud initiation factor homolog protein (19-211AA).

**Isotype:** IgG

**Form:** Liquid

**Purity:** > 95%

**Purification:** Purified by Protein G.

**Storage:** Aliquot and store at -20°C. Avoid exposure to light. Avoid repeated freeze/thaw cycles.

**UniProt Primary AC:** Q8TAL6 ([UniProt](#), [ExpAS](#))  
y

**Gene Symbol:** FIBIN

**GeneID:** [387758](#)

**OMIM:** [61785](#)

**KEGG:** hsa:387758

**String:** [9606.ENSP00000321962](#)

# Datasheet

Revision date: 08 Oct 2024



**BIOTECH ASSAY Inc.**

**Buffer:** 0.01 M PBS, pH 7.4, 0.03% Proclin-300 and 50% Glycerol.

**Note:** This product is for research use only.

BIOTECH ASSAY