

# Foot and Mouth Disease Virus Type O Antibody (Anti-FMDV-O) Rapid Test Kit

Catalogue No.:BTA11114

Foot and Mouth Disease Virus Type O Antibody Rapid Test Kit is a qualitative colloidal gold-based lateral flow immunochromatographic assay for the detection of Foot and Mouth Disease Virus Type O Antibody.

**Target:** Foot and Mouth Disease Virus Type O Antibody (Anti-FMDV-O)

**Reactivity:** Pig, Cow, Goat, Sheep

**Tested Applications:**LF

**Storage:** Store at 2-30 °C.

**Validity:** The validity for this kit is 24 months.

**Detection Method:** Immunochromatographic

**Assay Type:** Sandwich

**Assay Data:** Qualitative

**Sample Type:** Serum, plasma, whole blood.

**Note:** This product is for research use only.

Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different