

Human Hantavirus IgG/IgM (HV IgG/IgM) Rapid Test Kit

Catalogue No.:BTA11106

Hantavirus IgG/IgM (HV IgG/IgM) Rapid Test Kit is a qualitative colloidal gold-based lateral flow immunochromatographic assay for the detection of Hantavirus IgG/IgM (HV IgG/IgM).

Target: Hantavirus IgG/IgM (HV IgG/IgM)

Reactivity: Human

Tested Applications:LF

Storage: Store at 4°C to 30°C upon receipt.

Validity: The validity for this kit is 12 months.

Detection Method: Immunochromatographic

Assay Type: Capture

Assay Data: Qualitative

Sample Type: Serum

Note: This product is for research use only.

Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.