

## Human Treponema Pallidum (TP) Antibody Rapid Test Kit

Catalogue No.:BTA11109

Treponema Pallidum (TP) Antibody Rapid Test Kit is a qualitative colloidal gold-based lateral flow immunochromatographic assay for the detection of Treponema Pallidum (TP) Antibody.

**Target:** Treponema Pallidum (TP) Antibody

**Reactivity:** Human

**Tested Applications:**LF

**Storage:** Store at 4°C to 30°C upon receipt.

**Validity:** The validity for this kit is 12 months.

**Detection Method:** Immunochromatographic

**Assay Type:** Sandwich

**Assay Data:** Qualitative

**Note:** This product is for research use only.

Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.