

Melamine Rapid Test Kit

Catalogue No.:BTA11100

Melamine Rapid Test Kit is a qualitative colloidal gold-based lateral flow immunochromatographic assay for the detection of Melamine.

Target: Melamine

Reactivity: General

Tested Applications:LF

Storage: Store at 2-30 °C in a dry environment.

Sensitivity: 500 ppb (ng/ml)

Assay Type: Competitive Inhibition

Assay Data: Qualitative

Sample Type: Raw milk and milk.

Note: This product is for research use only.
Please note that our kits are optimised for detection of native samples, rather than recombinant proteins. We are unable to guarantee detection of recombinant proteins, as they may have different sequences or tertiary structures to the native protein.