

Mouse IgA Isotype Control Antibody

Catalogue No.:BTA14501

Mouse IgA Isotype Control Antibody is a negative control for measuring non-specific binding of Mouse IgA antibodies.

Target:	Mouse IgA
Clonality:	Monoclonal
Clone:	U332
Tested Applications:	FCM
Host:	Mouse
Conjugation:	Unconjugated
Isotype:	IgA
Form:	Liquid
Purity:	> 95%
Purification:	Purified by Protein A/G affinity chromatography.
Storage:	Store at 4 °C if the entire vial will be used within 1 week. Aliquot and store at -20 °C or -80 °C for long-term storage. Avoid repeated freeze-thaw cycles.
Buffer:	0.01 M PBS, pH 7.4.
Concentration:	1 mg/ml
Note:	This product is for research use only.