

Mouse IgG1 (FITC) / Mouse IgG1 (PE) Isotype Control

Catalogue No.:BTA14496

Mouse IgG1 (FITC) / Mouse IgG1 (PE) Isotype Control is a combination negative control cocktail conjugated to FITC/PE for measuring non-specific binding of Mouse IgG1 antibodies.

Target: Mouse IgG1

Clonality: Monoclonal

Tested Applications: FCM

Host: Mouse

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Conjugation: FITC, PE

Form: Lyophilized

Storage: Store lyophilized at 4 °C. After reconstitution, store undiluted at 4 °C. Do not freeze. This product is photosensitive and should be protected from light. If a precipitate is observed, we recommend microcentrifugation before use.

Buffer: Prior to lyophilization: PBS

Note: This product is for research use only.