

## Mouse IgG1 (FITC) / Mouse IgG2a (PE) Isotype Control

Catalogue No.:BTA14497

Mouse IgG1 (FITC) / Mouse IgG2a (PE) Isotype Control is a combination negative control cocktail conjugated to FITC/PE for measuring non-specific binding of Mouse IgG1 / Mouse IgG2a antibodies.

**Target:** Mouse IgG1/Mouse IgG2a

**Clonality:** Monoclonal

**Tested Applications:** FCM

**Host:** Mouse

**Recommended dilutions:** Optimal dilutions/concentrations should be determined by the end user.

**Conjugation:** FITC, PE

**Form:** Lyophilized

**Storage:** Store lyophilized at 4 °C. After reconstitution, store undiluted at 4 °C. Do not freeze. This product is photosensitive and should be protected from light. If a precipitate is observed, we recommend microcentrifugation before use.

**Buffer:** Prior to lyophilization: PBS

**Note:** This product is for research use only.