


Datasheet

Revision date: 07 Mar
2024

Mouse IgG1 Isotype Control (APC)

Catalogue No.:BTA14485



Fluorescence emission spectra of APC.

Mouse IgG1 Isotype Control is an Isotype Control conjugated to APC.

Target: Mouse IgG1

Clonality: Monoclonal

Clone: R828

Tested Applications: FCM

Host: Mouse

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Conjugation: APC

Excitation/Emission: 651/660

Laser Line: 647

Isotype: IgG, Kappa

Storage: Store in the dark at 2-8°C. Avoid exposure to light. Do not freeze.

Buffer: Stabilizing PBS solution containing 15 mM sodium azide.

Specificity: This mouse IgG1 kappa monoclonal antibody has unknown specificity and was chosen as an isotype control after screening on variety of resting, activated, live and fixed rat and human tissues.

Concentration: 0.1 mg/ml

Note: This product is for research use only.