

Mouse IgG1 Isotype Control

Catalogue No.:BTA14484

Mouse IgG1 Isotype Control is an Isotype Control.

Target:	Mouse IgG1
Clonality:	Monoclonal
Clone:	R828
Tested Applications:	ELISA, WB, IHC, IF/ICC, FCM, IP, FUNC
Host:	Mouse
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Conjugation:	Unconjugated
Isotype:	IgG, Kappa
Purity:	> 95% (SDS-PAGE)
Purification:	Purified by Protein A affinity chromatography.
Storage:	Store at 2-8°C. Do not freeze.
Buffer:	LE/AF PBS solution, 0.2 µm filter sterilized.
Endotoxin Level:	Endotoxin level is less than 0.01 EU/µg of the protein, as determined by the LAL test.
Specificity:	This mouse IgG1 kappa monoclonal antibody has unknown specificity and was chosen as an isotype control after screening on variety of resting, activated, live and fixed rat and human tissues.
Concentration:	1 mg/ml
Note:	This product is for research use only.