

Datasheet

Revision date: 09 Oct
2024

Mouse IgG2a Isotype Control

Catalogue No.:BTA14488

Mouse IgG2a Isotype Control is an Isotype Control.

Target: Mouse IgG2a

Clonality: Monoclonal

Clone: L584

Tested Applications: FCM

Host: Mouse

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Conjugation: Unconjugated

Isotype: IgG_{2a}

Purity: > 95% (SDS-PAGE)

Purification: Purified by Protein A affinity chromatography.

Storage: Store at 2-8°C. Do not freeze.

Buffer: PBS solution with 15 mM sodium azide.

Specificity: This mouse IgG2a monoclonal antibody reacts with undefined epitope on a plant pathogen.

Concentration: 1 mg/ml

Note: This product is for research use only.