


Datasheet

Revision date: 09 Oct
2024

Mouse IgG2b Isotype Control (FITC)

Catalogue No.:BTA14493



Fluorescence emission spectra of FITC.

Mouse IgG2b Isotype Control is an Isotype Control conjugated to FITC.

Target: Mouse IgG2b

Clonality: Monoclonal

Clone: C857

Tested Applications: FCM

Host: Mouse

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Conjugation: FITC

Excitation/Emission: 499/515

Laser Line: 488

Immunogen: KLH-coupled trinitrophenol.

Isotype: IgG_{2b} Kappa

Storage: Store in the dark at 2-8°C. Avoid exposure to light. Do not freeze.

Buffer: PBS solution with 15 mM sodium azide.

Specificity: This mouse IgG2b (kappa) monoclonal antibody reacts with an epitope irrelevant for a variety of resting, activated, live, and fixed human, mouse, and rat tissues.

Concentration: 1 mg/ml

Datasheet

Revision date: 09 Oct
2024



BIOTECH ASSAY Inc.

Note: This product is for research use only.

BIOTECH ASSAY