

Mouse IgG2b Isotype Control

Catalogue No.:BTA14490

Mouse IgG2b Isotype Control is an Isotype Control.

Target:	Mouse IgG2b
Clonality:	Monoclonal
Clone:	C857
Tested Applications:	ELISA, WB, IHC, IF/ICC, FCM, IP
Host:	Mouse
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Conjugation:	Unconjugated
Immunogen:	KLH-coupled trinitrophenol.
Isotype:	IgG _{2b} Kappa
Purity:	> 95% (SDS-PAGE)
Purification:	Purified by Protein A affinity chromatography.
Storage:	Store at 2-8°C. Do not freeze.
Buffer:	PBS solution with 15 mM sodium azide.
Specificity:	This mouse IgG2b (kappa) monoclonal antibody reacts with an epitope irrelevant for a variety of resting, activated, live, and fixed human, mouse, and rat tissues.
Concentration:	1 mg/ml
Note:	This product is for research use only.