

Rabbit Anti-Goat IgG Fab fragment Antibody (HRP)

Catalogue No.:BTA11022

Rabbit Anti-IgG Fab fragment Antibody (HRP) is a Rabbit Polyclonal antibody against IgG Fab fragment conjugated to HRP.

Target: Goat IgG Fab fragment

Clonality: Polyclonal

Reactivity: Goat

Tested Applications: ELISA

Host: Rabbit

Recommended dilutions: Optimal dilutions/concentrations should be determined by the end user.

Conjugation: HRP

Isotype: IgG

Form: Liquid

Purification: Purified using the caprylic acid ammonium sulfate (CA-AS) precipitation method.

Storage: Aliquot and store at -20°C. Avoid exposure to light. Avoid repeated freeze/thaw cycles.

Buffer: 0.01 M PBS, pH 7.4, 0.03% Proclin-300 and 50% glycerol.

Note: This product is for research use only.