

Rabbit Anti-IgG L Antibody (Biotin)

Catalogue No.:BTA11010

Rabbit Anti-IgG L Antibody (Biotin) is a Rabbit Polyclonal antibody against IgG L conjugated to Biotin.

Target:	IgG L
Clonality:	Polyclonal
Reactivity:	Human
Tested Applications:	ELISA
Host:	Rabbit
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Conjugation:	Biotin
Immunogen:	Human igg light chain protein (Native Protein).
Isotype:	IgG
Form:	Liquid
Purity:	> 95%
Purification:	Purified by Protein G.
Storage:	Aliquot and store at -20°C. Avoid exposure to light. Avoid repeated freeze/thaw cycles.
Buffer:	0.01 M PBS, pH 7.4, 0.03% Proclin-300 and 50% Glycerol.
Note:	This product is for research use only.