

Rabbit Anti-Mouse IgG Fc fragment Antibody (HRP)

Catalogue No.:BTA11025

Rabbit Anti-IgG Fc fragment Antibody (HRP) is a Rabbit Polyclonal antibody against IgG Fc fragment conjugated to HRP.

Target:	Mouse IgG Fc fragment
Clonality:	Polyclonal
Reactivity:	Mouse
Tested Applications:	ELISA
Host:	Rabbit
Recommended dilutions:	Optimal dilutions/concentrations should be determined by the end user.
Conjugation:	HRP
Immunogen:	Mouse IgG Fc fragment.
Isotype:	IgG
Form:	Liquid
Purification:	Purified using the caprylic acid ammonium sulfate (CA-AS) precipitation method.
Storage:	Store at -20°C or -80°C. Avoid exposure to light. Avoid repeated freeze/thaw cycles.
Buffer:	0.01 M PBS, pH 7.4, 0.03% Proclin-300 and 50% glycerol.
Note:	This product is for research use only.