

Datasheet

Revision date: 09 Oct
2024

Ubamatamab

Catalogue No.:BTA14515

This product is for research use only (RUO) by laboratories. It is not suitable for use in humans or as a medicine. Ubamatamab is a Chimeric Humanized Monoclonal antibody targeting Mucin-16, expressed in Mammalian cells.

Target:	Mucin-16
Clonality:	Monoclonal
Reactivity:	Human
Host:	CHO cells
Isotype:	IgG4 Kappa
Form:	Liquid
Purity:	> 95%
Purification:	Protein A or Protein G affinity purified.
Storage:	Store at 4°C for up to 2 weeks. Aliquot and store at -20°C or -80°C for up to 12 months.
UniProt Primary AC:	Q8WXI7 (UniProt , ExPAS) P07766 (UniProt , ExPAS)
Buffer:	PBS, pH 7.4
Concentration:	1 mg/ml
Note:	This product is for research use only. Not for human consumption, cosmetic, therapeutic or diagnostic use.