

## Xenopsin-Related Peptide 2

Catalogue No.:BTA11168

Xenopsin-Related Peptide 2 is a peptide.

<b>Target:</b>	Xenopsin-Related Peptide 2
<b>Conjugation:</b>	Unconjugated
<b>Form:</b>	Lyophilized
<b>Purity:</b>	> 95% (HPLC)
<b>Storage:</b>	Store dry at -20 °C. Avoid repeated freeze/thaw cycles.
<b>Molecular Weight:</b>	1193.47 Da
<b>Molecular Formula:</b>	$C_{60}H_{88}N_{16}O_{10}$
<b>Sequence:</b>	Phe-His-Pro-Lys-Arg-Pro-Trp-Ile-Leu
<b>Note:</b>	This product is for research use only. Not for human consumption, cosmetic, therapeutic or diagnostic use.